

- 
- 10ml syringe
    - Comes with 10 service kits
  - Must be used with the deep scrub swabs due to the extra strength gel (not included)
  - Apply the deep scrub swabs before the first and second (if you have enough) 15 minute interval for optimal results
    - 24K Gold Gel has bacteriostatic properties that help control oil bacteria
  - 24K Gold Gel is gentle for teeth and helps with sensitivity when used with the deep scrub swabs
  - Most luxurious Teeth Whitening Gel on the market
  - Contains approximately 0.06 Carats of 24K Gold Light Activated Nanoparticles to boost whitening results
    - Blended gel equivalent to 20% HP
  - 100% made in the US with Canadian gels
  - Contains the highest dental quality ingredients
  - Perfect for people with sensitive teeth, red wine, tea, coffee drinkers and smokers
    - Registered with Health Canada
    - Safe for veneers and crowns
  - Used and approved by dentists to efficiently and safely whiten smiles and remove stains from teeth
    - Keep refrigerated due to 24K Gold particles in gel
- Recommended Retail Pricing : \$140+